
as human, for human
PARAMOUNT BED

PARAMOUNT BED CO., LTD.
2-14-5 Higashisuna Koto-Ku, Tokyo 136-8670, Japan
Tel: +81 3 3648 2961 Fax: +81 3 3648 2420
http://www.paramount.co.jp/english/

Designs and specifications are subject to change without notice. All images of the product, its components and accessories may not reflect the final product. The images may reflect optional equipment and are subject to change without notice. All information presented herein is based on current data or design goals available at the time of posting and is subject to change without notice. Some equipment and/or features may be unavailable. Actual results of performance and other specifications may differ or vary with production models and may depend on selected models and/or accessories.

L152-1111-EN01

as human, for human
PARAMOUNT BED






Serie PA-90000

Cama Hospitalaria VIP

recuperación segura para su tranquiliad




PA-90000 VIP

Nuestra meta es proveer productos que ayudarán a la recuperación del paciente.
La cama PA-90000 VIP presenta un excelente diseño que mejorará la comodidad y tranquilidad del paciente. Damos un paso más al crear muebles que combinan visualmente con la cama, creando así un entorno positivo y propicio para el paciente.

Baja Altura de 290mm

Es importante mantener una posición firme en el suelo cuando los pacientes se levantan de la cama.
A una altura baja de 290 mm, incluso los usuarios con un físico pequeño pueden sentarse en el borde de la cama con los talones firmemente apoyados en el suelo..



Cama común
*Colchón común

Serie PA-90000 cama de altura baja
*Altura del colchón: 9cm

	Electric Bed Type
Código del producto	PA-93485V
Largo total	213.4 cm
Ancho total	116 cm
Altura total	76.1 to 114.1 cm
Altura base del colchón*	29 to 67 cm
Mattress Size	191 x 91 cm
Safe Working Load**	204kgf
Rueda	12,5 cm diá. ruedas dobles con sistema de bloqueo central
Ángulo de respaldo	0° to 70°
Ángulo de rodilla	0° to 25°

CE * Altura del suelo a la base del colchón.
**Carga máxima de operación (considerando el peso del paciente, accesorios y el colchón.

Accesorios adicionales

- Portasueros (KC-56A/B)
- Marco de tracción (KA-83)
- Soporte para tanque de oxígeno* (PC-5038)

*solo en la cabecera; puede contener 1,000L

DOUBLE-LOCKED SIDE-RAIL (PS-023V)



Optional Items

Paramount Bed offers a comprehensive variety of compatible furnishings to create a consistent and visually-appealing hospital ward environment.

- Wardrobe and Drawers





- Bedside Table (PF-3120AV)



- Sofa Bed (KF-Z09-01 series)



The backrest of the sofa bed can be reclined into a flat position.